




Step by step tick removal



1.

Locate the tick – they are usually found in the least hairy areas, particularly the face, neck, underbelly and inside of legs.

Choose the right tick twister based on the size of the tick.



2.

Put on gloves to avoid direct contact.

Approach the tick from the side using the tick twister.

Hook the ticks body.




3.

Lift the hook lightly whilst turning the twister.

After 2-3 rotations of the twister, the tick should release itself from your pet.

Put the tick in a container filled with alcohol and keep aside in case your vet needs to see this at a later date.



4.

Clean your pet's skin with antiseptic.

Thoroughly wash your hands.



Keep an eye on the area of your pet's skin where the tick was removed.

If the skin becomes irritated or infected, book an appointment to see a vet immediately.